

STERISHEET™ SMS Solutions



STERILIZATION WRAPS FOR HOSPITAL CSSDs
AND OUTSOURCED STERILIZATION UNITS

| Interleaved solutions |
| Single sheets | Bonded range |




STERISHEET™
by STERIMED



STERISHEET™

by STERIMED




STERISHEET – wide range of products for sterilization and sterility maintenance in patient care facilities.

Products you can trust, dedicated to prevent nosocomial infections.

STERISHEET PRODUCT RANGE


STERISHEET CREPE

100% Cellulose material




STERISHEET REINFORCED

Cellulose reinforced with synthetic binders




STERISHEET WET-LAID NON WOVEN

Combination of cellulose, synthetic fibers and binders



STERISHEET SMS NON WOVEN

100% synthetic fibers of polypropylene



- Can be delivered as:
- Single sheets
 - Interleaved, colour coded solution
 - for optimal safety in use
 - Sealed/bonded wraps
 - for quicker folding




STERISHEET SMS LINE

Our aim is to constantly develop our product range, improve supply chain, testing methods and monitor every need of our customers. Only continuous improvement and keeping up with latest trends bring us among the best companies in the world. Our products' range fits for any type of hospital sterilization method, size or weight of sterilized products.


Wide range of SMS products are now available to provide solutions for all hospital practices and regulations.



1 INTERLEAVED COMBINATIONS



2 SINGLE SHEETS



3 BONDED PRODUCTS

Use STERISHEET as interleaved, single sheets or bonded combination for the best results of your daily work in CSSD.

Standard sizes from 60x60 cm to 137x137 cm according to grade and grammage of the product. Other sizes available upon request.

- Excellent bacterial barrier
- Mechanical properties to withstand difficult conditions
- Memory free
- Linting free

STERISHEET SMS LINE










1 INTERLEAVED SOLUTIONS


According to recent hospital studies, double/sequential wrapping is recommended to improve sterility maintenance and to avoid particle contamination at opening. Therefore, interleaved products offer a better bacterial filtration efficiency and germ proofness.

What are the objectives of sequential wrapping?

- Create the most efficient tortuous path with the assembly of the different folds to avoid recontamination until opening
- Benefit from colour coding system
- Enable aseptic opening at the point of use

CHOOSE STERISHEET INTERLEAVED SOLUTION ACCORDING TO THE PACKS WEIGHT, PACKS LIFE CYCLE OR INTERNAL STANDARDS.

PACK WEIGHT	PRODUCT NAME	COLOUR
LOW 	Interleaved 340 SMS B / 340 SMS G	
	Interleaved 160 Soft G / 340 SMS B	
MEDIUM 	Interleaved 355 SMS B / 355 SMS G	
	Interleaved 260 Reinforced B / 355 SMS G	
HEAVY 	Interleaved 357 NW G / 363 SMS B	
	Interleaved 355 SMS G / 370 SMS B	



STERISHEET SMS LINE

2 SINGLE SHEETS

Expanded STERISHEET SMS portfolio gives you more options for your sterilization packaging.

- Two layers of meltblown (SMMS) for optimal barrier & filtration properties
- Pigments used for the offered SMS range of products are in compliance with European Cosmetics regulation
- Compliant with MDD 93/42/EEC amended by 2007/47/EC
- Compliant with EN ISO 11607:1-2017/EN 868:2-2017

5 GRADES, 2 COLOURS

PACK WEIGHT	PRODUCT NAME	GSM	Standard products	Available for specific needs	STERILIZATION
			COLOURS		
LOW 	STERISHEET 340 SMS	40			
MEDIUM 	STERISHEET 347 SMS	47			
	STERISHEET 355 SMS	55			
HEAVY 	STERISHEET 363 SMS	63			
	STERISHEET 370 SMS	70			



BACTERIAL BARRIER



TESTED



2 COLOURS

Don't put your patients under unnecessary risk, use sterilization wraps you can trust.

3 BONDED PRODUCTS

What are STERISHEET Bonded products?

Sealed combination of SMS & cellulose base product or SMS & SMS.

Two separate layers sealed together:

- Two colours – green/blue for easy defects detection
- Save your time during wrapping
- Double wrapping with single move
- Effective and efficient bacterial barrier


● STERISHEET DUO™

→ STERISHEET DUO™ Combination of SMS and high quality cellulose based products


→ Superior bacterial filtration efficiency thanks to cellulose

→ Cellulose helps to prevent wet pack issues

Warning: STERISHEET DUO™ products are not suitable for hydrogen peroxide (H₂O₂) gas plasma sterilization process





TRAY WEIGHT	PRODUCT NAME	PRODUCT DESCRIPTION	COLOUR	STERILIZATION
LOW 	STERISHEET Duo 1 (Crepe / SMS)	STERISHEET 160 Soft G / 340 SMS B		
MEDIUM 	STERISHEET Duo 2 (Reinforced / SMS)	STERISHEET 260 Reinforced B / 355 SMS G		
HEAVY 	STERISHEET Duo 3 (Non Woven / SMS)	STERISHEET 357 NW G / 363 SMS B		




● STERISHEET BONDED™

- Combination of SMS & SMS
- Suitable for most common types of sterilization process
- Excellent mechanical resistance thanks to 2 layers of SMS



SMS



SMS

TRAY WEIGHT	PRODUCT NAME	PRODUCT DESCRIPTION	COLOUR	STERILIZATION
LOW 	STERISHEET Bonded 2 (SMS / SMS)	STERISHEET 340 SMS B / 340 SMS G		
MEDIUM 	Sterisheet Bonded 3 (SMS / SMS)	Sterisheet 340 SMS G / 355 SMS B		
HEAVY 	STERISHEET Bonded 4 (SMS / SMS)	STERISHEET 355 SMS B / 355 SMS G		
SUPERHEAVY 	STERISHEET Bonded 5 (SMS / SMS)	STERISHEET 355 SMS G / 370 SMS B		

STERIMED, YOUR INFECTION CONTROL PARTNER:

- Largest range of sterilization packaging substrates in the world
- Sales in more than 70 countries
- Packaging education programmes for medical device industry and hospital staff

COMPLIANCE

Our products are compliant with strict regulations and standards, ensuring maximum protection for patients and hospital staff.

Our products conform with:

- EN ISO 11607-1:2017
- EN868-2:2017

We follow:

- Risk management
- Quality control
- ISO 13485 standards
- ISO 50001



PATIENT SAFETY FORUM EDUCATIONAL SEMINAR

We help professionals from hospitals & sterilization centres to share knowledge and get better understanding of material used to sterilize medical instruments in CSSDs.

www.patient-safety-forum.com



NEW DOWNLOAD STERIMED APP

A new platform for education, products selection and communication with STERIMED experts. Direct access to:

- Instructions For Use
- Technical Data
- Product Properties
- Frequently Asked Questions

Take an advantage of online tool for products' replenishment, samples ordering and trouble shooting.



STERISHEET™
by STERIMED

SALES OFFICES

AMERICAS

• **STERIMED Inc**
1301 Charleston Regional Parkway - Charleston 29492 - SC - South Carolina UNITED STATES
Phone: +1 843 388 80 80
Fax: +1 843 388 80 70

EUROPE, AFRICA & THE MIDDLE EAST

• **STERIMED SRO (Sterisheet™ product line)**
CT Park - Evropska 866 Modrice, 664 42 CZECH REPUBLIC
Phone: +420 547 424 934
Fax: +420 547 424 930

• **STERIMED SAS (Propypel™ & Ethypel™ product line)**
32 avenue Pierre Grenier 92100 Boulogne-Billancourt FRANCE
Phone: +33 157 75 92 75
Fax: +33 157 75 92 93

ASIA-PACIFIC

• **STERIMED (Suzhou) Co. Ltd**
Wujiang Economic & Technological Development Zone No. 1868, Guangming Road, Suzhou, Jiangsu, 215200, CHINA
Phone: +86 (0)512 6316 0176
Fax: +86 (0)512 6316 0151

CONTACT: contact@sterimed.fr

STERISHEET™ IS A BRAND OF



STERIMED MARKETING DEPARTMENT

32 avenue Pierre Grenier
92100 Boulogne-Billancourt - FRANCE
Phone: +33 157 75 92 75 - Fax: +33 1 57 75 92 93
contact@sterimed.fr



www.sterimed.fr



STERIMED™, ETHYPEL™, PROPYPEL™, STERISHEET™, ATH MEDICAL™, POLYBOND™, ETHYFLEX™, ETHYFORM™, CRYPTOPAK™, STERIWEB™, STERI FEUILLE™, 斯迪迈™, 舒卫安™, 泰洁安™, Стеримед™, Этипель™, ПРОПИПЕЛЬ™ are all brands of STERIMED.