


WRAPPING TECHNIQUES

> The packaging has 3 main functionalities:


- 1 • Provide the best bacterial barrier filtration
- 2 • Avoid re-contamination during storage and handling
- 3 • Facilitate aseptic opening

> The most common and recommended folding methods:

1 • Envelope folding



2 • Square folding



> Use double wrapping for the best results.

Double wrapping provides:

- 1 • The possibility to combine wraps of different composition and offer the best characteristics and benefits of each product component
- 2 • The possibility to have color coding, holes and tears are more easily detected with a color combination
- 3 • Enable aseptic opening at the point of use

TECHNICAL SHEET | STERISHEET WRAPING | Edited on september 2017.

STERISHEET®, A BRAND OF



STERIMED MARKETING DEPARTMENT
 32 avenue Pierre Grenier - 92100 Boulogne-Billancourt - FRANCE
 Phone: +33 157 75 92 75 - Fax: +33 1 57 75 92 93 - contact@sterimed.fr
www.sterimed.fr

