

Standard Basket Process Challenge Device (PCD)

- Stainless steel tray for testing of
- the ParaSure CEI
- Conforms to ISO 15883
- Easy to use –Easy to read
- Non toxic
- For use in Washer
- Disinfectors and Ultrasonic machines



Standard PCD (Washer/Disinfector Process Challenge Device)

- Container system for the ParaSure CEI (Cleaning Efficacy Indicator) creating a realistic and repeatable challenge. Process Challenge Device (PCD) for standard/general non-lumen instrumentation.

Order Details

PS1610

Standard Basket Process Challenge Device (PCD)





Registered Address:

ParaSure Limited Units B10/11 Oakendene Estate Bolney Rd Cowfold West Sussex RH13 8AZ England



Standard Basket Process Challenge Device (PCD)

TECHNICAL SPECIFICATIONS		
Name:	Standard Basket Process Challenge Device (PCD)	
Classification:	ISO 15883	
Storage conditions:	n/a	
RH:	n/a	
Lot no:	n/a	
Prod date:	n/a	
Exp date:	n/a	
Correct symbols:	n/a	
Indicator:	In conjunction with ParaSure CEI	
Disp/single use:	reusable	
Extra information: Non toxic	yes	

Assessment - Aid

Result	Assessment	Action needed	
	Reference indicator: Spot of test soil completely visible		
	INSUFFICIENT CLEANING PROCESS:		
	Spot of test soil completely visible	Do not release batch! Check load, temperature, dosage and time. Repeat test. In case of repeated insufficient cleaning process: Call service technician.	
0	Spot of test soil visible		
	Spot of test soil visible		
1	Spot of test soil slightly visible		
	SUCCESSFUL CLEANING PROCESS:		
A TOPY	Spot of test soil completely removed	Release batch. No further action needed.	

For further information about any of our products, please visit: www.parasure.co.uk





ParaSure Limited Units B10/11 Oakendene Estate Bolney Rd Cowfold West Sussex RH13 8AZ England

