



Manepa Hygienics, 628B Jordanstown Avenue,
Greenogue Business Park, Rathcoole, Co. Dublin
Tel: 01 4677600 Fax: 01 4105660
Website: www.manepa.com Email: support@manepa.com

TECHNICAL SPECIFICATION

ALD 29 – Concentrated Bathroom Cleaner

Highly fragranced, deep blue liquid developed for use around washroom areas. This multi-purpose detergent can be used on sinks, floors, and toilets, leaving surfaces clean and sparkling with a fresh and fragrant smell. Product is based on a blend of acid, surfactants, solvents, fragrance and water-soluble dyes.

Features include:

- Economical dilution rates.
- Deodorises and degreases in one.
- Removes body fats and limescale.
- Leaves room smelling clean and fresh.
- Removes the most stubborn soils found in washrooms such as soap scum.
- Kills odour producing bacteria eliminating the source of malodours.
- Not suitable for use on chrome/brass plated fittings unless thoroughly rinsed with clean water.
- All surfactants used are biodegradable and phosphate free.

PROPERTIES:

Appearance @ 20°C
pH (neat)
Odour
Specific Gravity @ 20°C
Flammability

SPECIFICATION:

Dark blue, clear liquid
>1
Fresh, floral fragrance
1.012 +/-0.005
non-flammable

DIRECTIONS FOR USE:

- * Use through appropriate dispenser.
- * Spray on surface to be cleaned, wipe off with a clean lint free cloth.
- * For heavy soil build up give product 2 minute contact time with surface before wiping.

STORAGE:

Store in original container and protect from extremes of temperature. Do not store in direct sunlight nor allow the product to freeze. Shelf life (in unopened original container) is a minimum of 18 months.

SAFETY ADVICE:

Under the new GHS Standard this product is classified as Irritant. Please refer to section 2.2 of the Safety Data Sheet for more information.

Manufactured to standard of BSEN 1276 testing. **PCS No: 100112**

PACK SIZE:

4x2L case (105 per pallet)



Manepa Hygienics, 628B Jordanstown Avenue,
Greenogue Business Park, Rathcoole, Co. Dublin
Tel: 01 4677600 Fax: 01 4105660
Website: www.manepa.com Email: support@manepa.com

TECHNICAL SPECIFICATION

ALD 29 – Concentrated Bathroom Cleaner

Reload No 9 is effective against the following bacteria / fungi / algae and yeasts at a dilution rate of 1:40 in 3 minutes:

BACTERIA

Salmonella typhosa
Salmonella enteritidis
Pseudomonas aeruginosa
Listeria monocytogenes
Shigella dysenteriae
Aerobacter aerogenes
Micrococcus pyrogenes
Streptococcus pyogenes
Streptococcus mutans
Clostridium tetani
Clostridium oedematiens
Klebsiella pneumoniae
Serratia marcescens
Candida albicans

Salmonella paratyphi
Salmonella dublin
Pseudomonas fluorescens
Desulphovibrio desulphuricans
Bacillus subtilis
Vibrio cholerae
Pneumococcus II
Streptococcus faecalis
Clostridium Welchii
Clostridium histolyticum
Klebsiella aerogenes
Enterobacter cloacae
Proteus vulgaris
Proteus mirabilis

FUNGI

Aspergillus niger
Myrothecium verucaria

Chaetomium globosum
Trichoderma viride

ALGAE

Chlorella vulgaris
Anabaena cylindrica

Stigeoclonium sp.
Oscillatoria tenuis

YEASTS

Saccharomyces cerevisiae
Endomycopsis albicans

Rhodotorula rubra

To eliminate MRSA Reload No.9 must be used at a dilution rate of 1:6.

ALD29 & BSEN 1276 TEST STANDARD

This British Standard test measures biocidal activity against *Escherichia coli*, amongst other organisms. The contact time is 5 minutes and the required kill is 5 log₁₀ reduction in organisms. Reload No.9 will pass this test at the standard dilution rate of 1:40.

Reload No.9 also passes this test in relation to other test organisms such as *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Proteus mirabilis*, *Streptococcus faecium* and *Saccharomyces cerevisiae*.